

CYLD Polyclonal Antibody

Cat No: HR1AP13004

For research use only

Overview

Product Name	CYLD Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	105kDa
GeneID?Human?	CYLD CYLD1 KIAA0849 HSPC057
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>CYLD lysine 63 deubiquitinase(CYLD) Homo sapiens This gene is encodes a cytoplasmic protein with three cytoskeletal-associated protein-glycine-conserved (CAP-GLY) domains that functions as a deubiquitinating enzyme. Mutations in this gene have been associated with cylindromatosis, multiple familial trichoepithelioma, and Brooke-Spiegler syndrome. Alternate transcriptional splice variants, encoding different isoforms, have been characterized. [provided by RefSeq, Jul 2008],</p>