

ST2A1 Polyclonal Antibody

Cat No: HR1AP12996

For research use only

Overview

Product Name	ST2A1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	31kDa
GeneID?Human?	SULT2A1 HST STD
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>sulfotransferase family 2A member 1(SULT2A1) Homo sapiens This gene encodes a member of the sulfotransferase family. Sulfotransferases aid in the metabolism of drugs and endogenous compounds by converting these substances into more hydrophilic water-soluble sulfate conjugates that can be easily excreted. This protein catalyzes the sulfation of steroids and bile acids in the liver and adrenal glands, and may have a role in the inherited adrenal androgen excess in women with polycystic ovary syndrome. [provided by RefSeq, Mar 2010],</p>