

# HHEX Polyclonal Antibody

Cat No: HR1AP12995

For research use only

## Overview

Product Name	HHEX Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	29kDa
GeneID?Human?	HHEX HEX PRH PRHX
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	hematopoietically expressed homeobox(HHEX) Homo sapiens This gene encodes a member of the homeobox family of transcription factors, many of which are involved in developmental processes. Expression in specific hematopoietic lineages suggests that this protein may play a role in hematopoietic differentiation. [provided by RefSeq, Jul 2008],