

MUC3A Polyclonal Antibody

Cat No: HR1AP12994

For research use only

Overview

Product Name	MUC3A Polyclonal Antibody
Source	Rabbit
Applications	IHC-p
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	279kDa
GeneID?Human?	MUC3A MUC3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>mucin 3A, cell surface associated(MUC3A) Homo sapiens The mucin genes encode epithelial glycoproteins, some of which are secreted and some membrane bound. Each of the genes contains at least one large domain of tandemly repeated sequence that encodes the peptide sequence rich in serine and/or threonine residues, which carries most of the O-linked glycosylation (Gendler and Spicer, 1995 [PubMed 7778880]).[supplied by OMIM, Aug 2008].</p>