

RAGP1 Polyclonal Antibody

Cat No: HR1AP12991

For research use only

Overview

Product Name	RAGP1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	64kDa
GeneID?Human?	RANGAP1 KIAA1835 SD
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	Ran GTPase activating protein 1(RANGAP1) Homo sapiens This gene encodes a protein that associates with the nuclear pore complex and participates in the regulation of nuclear transport. The encoded protein interacts with Ras-related nuclear protein 1 (RAN) and regulates guanosine triphosphate (GTP)-binding and exchange. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2013],