

EIF1 Polyclonal Antibody

Cat No: HR1AP12990

For research use only

Overview

Product Name	EIF1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	12kDa
GeneID?Human?	EIF1 SUI1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:Necessary for scanning and involved in initiation site selection. Promotes the assembly of 48S ribosomal complexes at the authentic initiation codon of a conventional capped mRNA.,similarity:Belongs to the SUI1 family.,subunit:Interacts with RENT2.,