

# ETS2 Polyclonal Antibody

Cat No: HR1AP12987

For research use only

## Overview

Product Name	ETS2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	51kDa
GeneID?Human?	ETS2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	ETS proto-oncogene 2, transcription factor(ETS2) Homo sapiens This gene encodes a transcription factor which regulates genes involved in development and apoptosis. The encoded protein is also a protooncogene and shown to be involved in regulation of telomerase. A pseudogene of this gene is located on the X chromosome. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2012],