

DIAP2 Polyclonal Antibody

Cat No: HR1AP12984

For research use only

Overview

Product Name	DIAP2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	121kDa
GeneID?Human?	DIAPH2 DIA
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>diaphanous related formin 2(DIAPH2) Homo sapiens The product of this gene belongs to the diaphanous subfamily of the formin homology family of proteins. This gene may play a role in the development and normal function of the ovaries. Defects in this gene have been linked to premature ovarian failure 2. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jul 2008],</p>