

PRKRA Polyclonal Antibody

Cat No: HR1AP12985

For research use only

Overview

Product Name	PRKRA Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	34kDa
GeneID?Human?	PRKRA PACT RAX HSD-14 HSD14
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	protein activator of interferon induced protein kinase EIF2AK2(PRKRA) Homo sapiens This gene encodes a protein kinase activated by double-stranded RNA which mediates the effects of interferon in response to viral infection. Mutations in this gene have been associated with dystonia. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Nov 2008],