

GON2 Polyclonal Antibody

Cat No: HR1AP12983

For research use only

Overview

Product Name	GON2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	13kDa
GeneID?Human?	GNRH2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>gonadotropin releasing hormone 2(GNRH2) Homo sapiens This gene encodes a secreted peptide hormone and member of the gonadotropin-releasing hormone (GnRH) family of proteins. The encoded protein regulates reproductive function by stimulating the production and release of the gonadotropins follicle-stimulating hormone (FSH) and luteinizing hormone (LH). The encoded protein may inhibit endometrial, ovarian, prostate, and breast cancer cell proliferation. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2015].</p>