

S5A1 Polyclonal Antibody

Cat No: HR1AP12976

For research use only

Overview

Product Name	S5A1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	28kDa
GeneID?Human?	SRD5A1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	steroid 5 alpha-reductase 1(SRD5A1) Homo sapiens Steroid 5-alpha-reductase (EC 1.3.99.5) catalyzes the conversion of testosterone into the more potent androgen, dihydrotestosterone (DHT). Also see SRD5A2 (MIM 607306).[supplied by OMIM, Mar 2008],