

SRTD1 Polyclonal Antibody

Cat No: HR1AP12971

For research use only

Overview

| | |
|-----------------------|--|
| Product Name | SRTD1 Polyclonal Antibody |
| Source | Rabbit |
| Applications | WB,ELISA |
| Species Reactivity | Human |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 25kDa |
| GeneID?Human? | SERTAD1 SE11 |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | |
| Background | function:Acts at E2F-responsive promoters to integrate signals provided by PHD- and/or bromodomain-containing transcription factors. Stimulates E2F-1/DP-1 transcriptional activity. Renders the activity of cyclin D1/CDK4 resistant to the inhibitory effects of p16(INK4a).,similarity:Contains 1 SERTA domain.,subunit:Interacts with the PHD-bromodomain of TIF1, TRIM28/TIF1B and p300/CBP. Binds to DP1. Also interacts with CDK4., |