

RPRM Polyclonal Antibody

Cat No: HR1AP12968

For research use only

Overview

Product Name	RPRM Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	11kDa
GeneID?Human?	RPRM
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:May be involved in the regulation of p53-dependent G2 arrest of the cell cycle. Seems to induce cell cycle arrest by inhibiting CDC2 activity and nuclear translocation of the CDC2 cyclin B1 complex.,induction:By TP53/p53, following X-ray irradiation.,miscellaneous:'Reprimo' signifies stop/repress.,similarity:Belongs to the reprimo family.,