



PODO Polyclonal Antibody

Cat No: HR1AP12966

For research use only

Overview

Product Name	PODO Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	42kDa
GeneID?Human?	NPHS2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	NPHS2, podocin(NPHS2) Homo sapiens This gene encodes a protein that plays a role in the regulation of glomerular permeability. Mutations in this gene cause steroid-resistant nephrotic syndrome. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2014],