

PCY1A Polyclonal Antibody

Cat No: HR1AP12965

For research use only

Overview

Product Name	PCY1A Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	40kDa
GeneID?Human?	PCYT1A CTPCT PCYT1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	phosphate cytidyltransferase 1, choline, alpha(PCYT1A) Homo sapiens This gene belongs to the cytidyltransferase family and is involved in the regulation of phosphatidylcholine biosynthesis. Mutations in this gene are associated with spondylometaphyseal dysplasia with cone-rod dystrophy. Alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Aug 2015],