

MGRN1 Polyclonal Antibody

Cat No: HR1AP12962

For research use only

Overview

Product Name	MGRN1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	60kDa
GeneID?Human?	MGRN1 KIAA0544 RNF156
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	mahogunin ring finger 1(MGRN1) Homo sapiens Mahogunin (MGRN1) is a C3HC4 RING-containing protein with E3 ubiquitin ligase activity in vitro.[supplied by OMIM, Apr 2004],