

PROL4 Polyclonal Antibody

Cat No: HR1AP12960

For research use only

Overview

Product Name	PROL4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	14kDa
GeneID?Human?	PRR4 LPRP PROL4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>proline rich 4 (lacrimal)(PRR4) Homo sapiens This gene encodes a member of the proline-rich protein family that lacks a conserved repetitive domain. This protein may play a role in protective functions in the eye. Alternative splicing result in multiple transcript variants. Read-through transcription also exists between this gene and the upstream PRH1 (proline-rich protein HaeIII subfamily 1) gene. [provided by RefSeq, Feb 2011],</p>