

IGKC Polyclonal Antibody

Cat No: HR1AP12957

For research use only

Overview

Product Name	IGKC Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	11kDa
GeneID?Human?	IGKC
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	miscellaneous:The EU sequence has the INV (3) allotypic marker, Ala-45 and Val-83. The ROY sequence has the INV (1,2) allotypic marker, Ala-45 and Leu-83.,similarity:Contains 1 Ig-like (immunoglobulin-like) domain.,