

## IMMT Polyclonal Antibody

Cat No: HR1AP12956

For research use only

### Overview

Product Name	IMMT Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	83kDa
GeneID?Human?	IMMT HMP MINOS2 PIG4 PIG52
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	Component of the MICOS complex, a large protein complex of the mitochondrial inner membrane that plays crucial roles in the maintenance of crista junctions, inner membrane architecture, and formation of contact sites to the outer membrane. Plays an important role in the maintenance of the MICOS complex stability and the mitochondrial cristae morphology (PubMed: 22114354, PubMed: 25781180).