

## CGT Polyclonal Antibody

Cat No: HR1AP12942

For research use only

### Overview

Product Name	CGT Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	59kDa
GeneID?Human?	UGT8 CGT UGT4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>UDP glycosyltransferase 8(UGT8) Homo sapiens The protein encoded by this gene belongs to the UDP-glycosyltransferase family. It catalyzes the transfer of galactose to ceramide, a key enzymatic step in the biosynthesis of galactocerebrosides, which are abundant sphingolipids of the myelin membrane of the central and peripheral nervous systems. Alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Sep 2011].</p>