

CEGT Polyclonal Antibody

Cat No: HR1AP12936

For research use only

Overview

Product Name	CEGT Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	43kDa
GeneID?Human?	UGCG
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>UDP-glucose ceramide glucosyltransferase(UGCG) Homo sapiens This gene encodes an enzyme that catalyzes the first glycosylation step in the biosynthesis of glycosphingolipids, which are membrane components containing lipid and sugar moieties. The product of this reaction is glucosylceramide, which is the core structure of many glycosphingolipids. [provided by RefSeq, Dec 2014],</p>