

TM9S1 Polyclonal Antibody

Cat No: HR1AP12933

For research use only

Overview

| | |
|-----------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name | TM9S1 Polyclonal Antibody |
| Source | Rabbit |
| Applications | WB,ELISA |
| Species Reactivity | Human,Mouse,Rat |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 66kDa |
| GeneID?Human? | TM9SF1 |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | |
| Background | similarity:Belongs to the nonaspanin (TM9SF) family.,tissue specificity:Expressed in lung, pancreas, kidney, liver, placenta, skeletal muscle, heart and brain. The amount in skeletal muscle, heart and brain were considerably lower than in the other tissues., |