

## RARR2 Polyclonal Antibody

Cat No: HR1AP12932

For research use only

### Overview

Product Name	RARR2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	17kDa
GeneID?Human?	RARRES2 TIG2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	retinoic acid receptor responder 2(RARRES2) Homo sapiens This gene encodes a secreted chemotactic protein that initiates chemotaxis via the ChemR23 G protein-coupled seven-transmembrane domain ligand. Expression of this gene is upregulated by the synthetic retinoid tazarotene and occurs in a wide variety of tissues. The active protein has several roles, including that as an adipokine and as an antimicrobial protein with activity against bacteria and fungi. [provided by RefSeq, Nov 2014],