

## SEM4C Polyclonal Antibody

Cat No: HR1AP12929

For research use only

### Overview

Product Name	SEM4C Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	91kDa
GeneID?Human?	SEMA4C KIAA1739 SEMAI UNQ5855/PRO34487
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	similarity:Belongs to the semaphorin family.,similarity:Contains 1 Ig-like C2-type (immunoglobulin-like) domain.,similarity:Contains 1 PSI domain.,similarity:Contains 1 Sema domain.,subunit:Interacts with GIPC PDZ domain. Interacts with NCDN.,