

P2RX2 Polyclonal Antibody

Cat No: HR1AP12924

For research use only

Overview

Product Name	P2RX2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	51kDa
GeneID?Human?	P2RX2 P2X2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>purinergic receptor P2X 2(P2RX2) Homo sapiens The product of this gene belongs to the family of purinoceptors for ATP. This receptor functions as a ligand-gated ion channel. Binding to ATP mediates synaptic transmission between neurons and from neurons to smooth muscle. Multiple transcript variants encoding distinct isoforms have been identified for this gene. [provided by RefSeq, Aug 2013],</p>