

## ESM1 Polyclonal Antibody

Cat No: HR1AP12917

For research use only

### Overview

Product Name	ESM1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	20kDa
GeneID?Human?	ESM1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>endothelial cell specific molecule 1(ESM1) Homo sapiens This gene encodes a secreted protein which is mainly expressed in the endothelial cells in human lung and kidney tissues. The expression of this gene is regulated by cytokines, suggesting that it may play a role in endothelium-dependent pathological disorders. The transcript contains multiple polyadenylation and mRNA instability signals. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2008],</p>