

MOAP1 Polyclonal Antibody

Cat No: HR1AP12914

For research use only

Overview

Product Name	MOAP1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	38kDa
GeneID?Human?	MOAP1 PNMA4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>modulator of apoptosis 1(MOAP1) Homo sapiens The protein encoded by this gene was identified by its interaction with apoptosis regulator BAX protein. This protein contains a Bcl-2 homology 3 (BH3)-like motif, which is required for the association with BAX. When overexpressed, this gene has been shown to mediate caspase-dependent apoptosis. [provided by RefSeq, Jul 2008],</p>