

THYN1 Polyclonal Antibody

Cat No: HR1AP12912

For research use only

Overview

Product Name	THYN1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	24kDa
GeneID?Human?	THYN1 THY28 HSPC144 MDS012 My0054
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	thymocyte nuclear protein 1(THYN1) Homo sapiens This gene encodes a protein that is highly conserved among vertebrates and plant species and may be involved in the induction of apoptosis. Alternatively spliced transcript variants encoding different isoforms have been described. [provided by RefSeq, Jul 2008],