

# SCG1 Polyclonal Antibody

Cat No: HR1AP12911

For research use only

## Overview

Product Name	SCG1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	74kDa
GeneID?Human?	CHGB SCG1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	chromogranin B(CHGB) Homo sapiens This gene encodes a tyrosine-sulfated secretory protein abundant in peptidergic endocrine cells and neurons. This protein may serve as a precursor for regulatory peptides. [provided by RefSeq, Jan 2009],