

TXLNA Polyclonal Antibody

Cat No: HR1AP12899

For research use only

Overview

Product Name	TXLNA Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	60kDa
GeneID?Human?	TXLNA TXLN
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>caution:Was originally (PubMed:8327514) thought to be a high molecular weight interleukin (IL-14 or IL14).,function:May be involved in intracellular vesicle traffic and potentially in calcium-dependent exocytosis in neuroendocrine cells.,similarity:Belongs to the taxilin family.,subunit:Binds to the C-terminal coiled coil region of syntaxin family members STX1A, STX3A and STX4A, but not when these proteins are complexed with SNAP25, VAMP2 or STXBP1, suggesting that it interacts with syntaxins that do not form the SNARE complex.,tissue specificity:Ubiquitous, with much higher expression in heart, kidney, liver and pancreas.,</p>