

ESYT1 Polyclonal Antibody

Cat No: HR1AP12895

For research use only

Overview

Product Name	ESYT1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	121kDa
GeneID?Human?	ESYT1 FAM62A KIAA0747 MBC2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:May play a role as calcium-regulated intrinsic membrane protein.,similarity:Belongs to the extended synaptotagmin family.,similarity:Contains 5 C2 domains.,subcellular location:Localizes to intracellular membranes.,tissue specificity:Widely expressed.,