

PKP3 Polyclonal Antibody

Cat No: HR1AP12896

For research use only

Overview

Product Name	PKP3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	87kDa
GeneID?Human?	PKP3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>plakophilin 3(PKP3) Homo sapiens This gene encodes a member of the arm-repeat (armadillo) and plakophilin gene families. Plakophilin proteins contain numerous armadillo repeats, localize to cell desmosomes and nuclei, and participate in linking cadherins to intermediate filaments in the cytoskeleton. This protein may act in cellular desmosome-dependent adhesion and signaling pathways. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Dec 2014],</p>