

AF1L2 Polyclonal Antibody

Cat No: HR1AP12893

For research use only

Overview

Product Name	AF1L2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	89kDa
GeneID?Human?	AFAP1L2 KIAA1914 XB130
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:May play a role in a signaling cascade by enhancing the kinase activity of SRC. Contributes to SRC-regulated transcription activation.,PTM:Tyrosine phosphorylated (by SRC).,similarity:Contains 2 PH domains.,subunit:Interacts with SRC. Interacts with LCK when tyrosine phosphorylated.,tissue specificity:Detected in spleen and thyroid, and at lower levels in kidney, brain, lung and pancreas.,