

CSKI2 Polyclonal Antibody

Cat No: HR1AP12894

For research use only

Overview

Product Name	CSKI2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	132kDa
GeneID?Human?	CASKIN2 KIAA1139
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>CASK interacting protein 2(CASKIN2) Homo sapiens This gene encodes a large protein that contains six ankyrin repeats, as well as a Src homology 3 (SH3) domain and two sterile alpha motif (SAM) domains, which may be involved in protein-protein interactions. The C-terminal portion of this protein is proline-rich and contains a conserved region. A related protein interacts with calcium/calmodulin-dependent serine protein kinase (CASK). Alternative splicing results in multiple transcript variants. [provided by RefSeq, Mar 2013],</p>