

LIMA1 Polyclonal Antibody

Cat No: HR1AP12892

For research use only

Overview

Product Name	LIMA1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	83kDa
GeneID?Human?	LIMA1 EPLIN SREBP3 PP624
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	LIM domain and actin binding 1 (LIMA1) Homo sapiens This gene encodes a cytoskeleton-associated protein that inhibits actin filament depolymerization and cross-links filaments in bundles. It is downregulated in some cancer cell lines. Alternatively spliced transcript variants encoding different isoforms have been described for this gene, and expression of some of the variants maybe independently regulated. [provided by RefSeq, Aug 2011],