

ARI3A Polyclonal Antibody

Cat No: HR1AP12890

For research use only

Overview

Product Name	ARI3A Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	65kDa
GeneID?Human?	ARID3A DRIL1 DRIL3 DRX E2FBP1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>AT-rich interaction domain 3A(ARID3A) Homo sapiens This gene encodes a member of the ARID (AT-rich interaction domain) family of DNA binding proteins. It was found by homology to the Drosophila dead ringer gene, which is important for normal embryogenesis. Other ARID family members have roles in embryonic patterning, cell lineage gene regulation, cell cycle control, transcriptional regulation, and possibly in chromatin structure modification. [provided by RefSeq, Jul 2008],</p>