

## SEM4B Polyclonal Antibody

Cat No: HR1AP12888

For research use only

### Overview

Product Name	SEM4B Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	91kDa
GeneID?Human?	SEMA4B KIAA1745 SEMAC UNQ749/PRO1480
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:Inhibits axonal extension by providing local signals to specify territories inaccessible for growing axons.,similarity:Belongs to the semaphorin family.,similarity:Contains 1 Ig-like C2-type (immunoglobulin-like) domain.,similarity:Contains 1 PSI domain.,similarity:Contains 1 Sema domain.,