

DC1L1 Polyclonal Antibody

Cat No: HR1AP12887

For research use only

Overview

Product Name	DC1L1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	57kDa
GeneID?Human?	DYNC1LI1 DNCLI1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>dynein cytoplasmic 1 light intermediate chain 1(DYNC1LI1) Homo sapiens The protein encoded by this gene belongs to light intermediate subunit family, whose members are components of the multiprotein cytoplasmic dynein complex, which is involved in intracellular trafficking and chromosome segregation during mitosis. The protein plays a role in moving the spindle assembly checkpoint (SAC) from kinetochores to spindle poles. The protein may also mediate binding to other cargo molecules to facilitate intracellular vesicle trafficking. [provided by RefSeq, Jul 2016],</p>