

MATR3 Polyclonal Antibody

Cat No: HR1AP12884

For research use only

Overview

Product Name	MATR3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	93kDa
GeneID?Human?	MATR3 KIAA0723
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>matrin 3(MATR3) Homo sapiens This gene encodes a nuclear matrix protein, which is proposed to stabilize certain messenger RNA species. Mutations of this gene are associated with distal myopathy 2, which often includes vocal cord and pharyngeal weakness. Alternatively spliced transcript variants, including read-through transcripts composed of the upstream small nucleolar RNA host gene 4 (non-protein coding) and matrin 3 gene sequence, have been identified. Pseudogenes of this gene are located on chromosomes 1 and X. [provided by RefSeq, Aug 2013],</p>