

ANS1A Polyclonal Antibody

Cat No: HR1AP12882

For research use only

Overview

Product Name	ANS1A Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	124kDa
GeneID?Human?	ANKS1A ANKS1 KIAA0229 ODIN
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:May play a negative role in growth factor receptor signaling pathways.,PTM:Phosphorylated on tyrosine residues in response to EGF and PDGF.,similarity:Contains 1 PID domain.,similarity:Contains 2 SAM (sterile alpha motif) domains.,similarity:Contains 6 ANK repeats.,subcellular location:Cytoplasmic before and after growth factor treatment.,tissue specificity:Widely expressed (at protein level),.