

DNMBP Polyclonal Antibody

Cat No: HR1AP12880

For research use only

Overview

Product Name	DNMBP Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	173kDa
GeneID?Human?	DNMBP KIAA1010
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>dynammin binding protein(DNMBP) Homo sapiens This gene encodes a protein belonging to the guanine nucleotide exchange factor family, and which regulates the configuration of cell junctions. It contains multiple binding sites for dynammin and thus links dynammin to actin regulatory proteins. Polymorphisms in this gene have been linked to Alzheimer's disease in some populations, though there are conflicting reports of such linkages in other populations. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2015],</p>