

DLG1 Polyclonal Antibody

Cat No: HR1AP12878

For research use only

Overview

Product Name	DLG1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	99kDa
GeneID?Human?	DLG1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	discs large MAGUK scaffold protein 1(DLG1) Homo sapiens This gene encodes a multi-domain scaffolding protein that is required for normal development. This protein may have a role in septate junction formation, signal transduction, cell proliferation, synaptogenesis and lymphocyte activation. Several alternatively spliced transcript variants encoding different isoforms have been described for this gene, but the full-length nature of some of the variants is not known. [provided by RefSeq, Feb 2011].