

## HNRPR Polyclonal Antibody

Cat No: HR1AP12876

For research use only

### Overview

Product Name	HNRPR Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	69kDa
GeneID?Human?	HNRNPR HNRPR
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	heterogeneous nuclear ribonucleoprotein R(HNRNPR) Homo sapiens This gene encodes an RNA-binding protein that is a member of the spliceosome C complex, which functions in pre-mRNA processing and transport. The encoded protein also promotes transcription at the c-fos gene. Alternative splicing results in multiple transcript variants. There are pseudogenes for this gene on chromosomes 4, 11, and 10. [provided by RefSeq, Jul 2014],