

KIRR1 Polyclonal Antibody

Cat No: HR1AP12873

For research use only

Overview

Product Name	KIRR1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	83kDa
GeneID?Human?	KIRREL KIRREL1 NEPH1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>kin of IRRE like (Drosophila)(KIRREL) Homo sapiens NEPH1 is a member of the nephrin-like protein family, which includes NEPH2 (MIM 607761) and NEPH3 (MIM 607762). The cytoplasmic domains of these proteins interact with the C terminus of podocin (NPHS2; MIM 604766), and the genes are expressed in kidney podocytes, cells involved in ensuring size- and charge-selective ultrafiltration (Sellin et al., 2003 [PubMed 12424224]).[supplied by OMIM, Mar 2008],</p>