

CDV3 Polyclonal Antibody

Cat No: HR1AP12866

For research use only

Overview

Product Name	CDV3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	28kDa
GeneID?Human?	CDV3 H41
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	similarity:Belongs to the CDV3 family.,tissue specificity:Expression levels correlate with those of HER-2/neu in breast cancer cells.,