

## CASL Polyclonal Antibody

Cat No: HR1AP12862

For research use only

### Overview

Product Name	CASL Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	91kDa
GeneID?Human?	NEDD9 CASL CASS2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>neural precursor cell expressed, developmentally down-regulated 9(NEDD9) Homo sapiens The protein encoded by this gene is a member of the CRK-associated substrates family. Members of this family are adhesion docking molecules that mediate protein-protein interactions for signal transduction pathways. This protein is a focal adhesion protein that acts as a scaffold to regulate signaling complexes important in cell attachment, migration and invasion as well as apoptosis and the cell cycle. This protein has also been reported to have a role in cancer metastasis. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2012],</p>