

# TKT Polyclonal Antibody

Cat No: HR1AP12861

For research use only

## Overview

Product Name	TKT Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	68kDa
GeneID?Human?	TKT
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	transketolase(TKT) Homo sapiens This gene encodes a thiamine-dependent enzyme which plays a role in the channeling of excess sugar phosphates to glycolysis in the pentose phosphate pathway. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Apr 2012],