

LSR Polyclonal Antibody

Cat No: HR1AP12858

For research use only

Overview

Product Name	LSR Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	71kDa
GeneID?Human?	LSR LISCH
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>caution:In contrast to the rodent orthologous protein, it is longer in N-terminus and no signal sequence is detected by any prediction method.,caution:The sequence shown here is derived from an Ensembl automatic analysis pipeline and should be considered as preliminary data.,function:Probable role in the clearance of triglyceride-rich lipoprotein from blood. Binds chylomicrons, LDL and VLDL in presence of free fatty acids and allows their subsequent uptake in the cells.,similarity:Belongs to the immunoglobulin superfamily. LISCH7 family.,similarity:Contains 1 Ig-like V-type (immunoglobulin-like) domain.,subunit:Homotrimer or homotetramer.,</p>