

CLH1 Polyclonal Antibody

Cat No: HR1AP12856

For research use only

Overview

Product Name	CLH1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	184kDa
GeneID?Human?	CLTC CLH17 CLTCL2 KIAA0034
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	clathrin heavy chain(CLTC) Homo sapiens Clathrin is a major protein component of the cytoplasmic face of intracellular organelles, called coated vesicles and coated pits. These specialized organelles are involved in the intracellular trafficking of receptors and endocytosis of a variety of macromolecules. The basic subunit of the clathrin coat is composed of three heavy chains and three light chains. [provided by RefSeq, Jul 2008],