

ITSN2 Polyclonal Antibody

Cat No: HR1AP12857

For research use only

Overview

Product Name	ITSN2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	186kDa
GeneID?Human?	ITSN2 KIAA1256 SH3D1B SWAP
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>intersectin 2(ITSN2) Homo sapiens This gene encodes a cytoplasmic protein which contains SH3 domains. This protein is a member of a family of proteins involved in clathrin-mediated endocytosis. Intersectin 2 is thought to regulate the formation of clathrin-coated vesicles and also may function in the induction of T cell antigen receptor (TCR) endocytosis. Alternatively spliced transcript variants have been found for this gene that encode three distinct isoforms. Additional variants have been found but their full length nature has not been determined. [provided by RefSeq, Jul 2008],</p>